

Victor Import Module

OVERVIEW

Abstract

The PLA 2.1 Import Module for Victor Plate Readers enables users to create biological assays in PLA 2.1 containing measurement values from these plate readers, provided as .txt files.

Features

- Fast, secure and reliable data transfer to PLA.
- Option to import several plates at once.
- All imported data is fully retraceable (audit trail).

TECHNICAL SPECIFICATION

Type

Import Module

System Requirements

- PLA 2.0 or PLA 2.1

Supported File Formats

- Victor .txt

Version

Version 2.1

(downgrade rights included)

Supported Systems

- Victor plate readers

Supported Languages

English, German

DELIVERY AND SUPPORT

Release Date

June 16, 2014

Scope of Delivery

- Software (CD or download)
- User manual (PDF)

Support and Upgrades

You have access to our PLA support platform at: support.bioassay.de

Product No.

23039

Lifecycle

We will provide support for PLA 2.1 and all related modules until May 8, 2023.

This module is covered by our PLA Support Contract. The first year is included and the contract can be renewed year by year. To find out more about its benefits, please go to: www.bioassay.de/support-contract

AVAILABLE MODULES FOR PLA 2.1

- Biacore™ Import Module
- Bio-Plex® Import Module
- Dynex Revelation™ Import Module
- Gen5™ Import Module
- BioTek® Instruments KC4™ Import Module
- Meso Scale Diagnostics® Import Module
- Microsoft Excel™ Import Module
- PerkinElmer® EnVision Import Module
- PerkinElmer® Victor Import Module
- SoftMax® Pro Import Module
- Tecan Magellan™ Import Module

Modules for other data sources will be developed on demand.