



Bio-Plex® Import Module

OVERVIEW

Abstract

The PLA 2.1 Import Module for Bio-Plex (Bio-Rad Laboratories, Inc) enables users to import XML-files created with Bio-Plex Manager 4.1 into PLA 2.1 for further analysis.

Features

- Fast, secure and reliable data transfer to PLA.
- All imported data is fully retraceable (audit trail).
- Option to create multiple assays from a single source file.

TECHNICAL SPECIFICATION

Type

Import Module

System Requirements

- PLA 2.0 or PLA 2.1

Supported File Formats

- .xml

Version

Version 1.4

(downgrade rights included)

Supported Systems

- Bio-Plex 100 and 200 systems
- Bio-Plex instruments

Supported Languages

English, German

DELIVERY AND SUPPORT

Release Date

March 13, 2015

Scope of Delivery

- Software (CD or download)
- User manual (PDF)

Support and Upgrades

You have access to our PLA support platform at: support.bioassay.de

Product No.

23036

Lifecycle

We will provide support for PLA 2.1 and all related modules until May 8, 2023.

This module is covered by our PLA Support Contract. The first year is included and the contract can be renewed year by year. To find out more about its benefits, please go to: www.bioassay.de/support-contract

AVAILABLE MODULES FOR PLA 2.1

- Biacore™ Import Module
- Bio-Plex® Import Module
- Dynex Revelation™ Import Module
- Gen5™ Import Module
- BioTek® Instruments KC4™ Import Module
- Meso Scale Diagnostics® Import Module
- Microsoft Excel™ Import Module
- PerkinElmer® EnVision Import Module
- PerkinElmer® Victor Import Module
- SoftMax® Pro Import Module
- Tecan Magellan™ Import Module

Modules for other data sources will be developed on demand.